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
### Motivation

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

## Heating Microscope

### Introduction

The heating microscope enables characteristic temperatures of materials to be measured, such as melting temperatures and sintering temperatures. Furthermore, it enables deformation processes to be observed as well as surface energy to be measured by determining the contact angle of the molten material.

### Measurement Principle

**Fig. 1:** Diagram of the heating microscope.


The apparatus consists of an oven, a light source and a microscope with an attached CCD-camera.

Small samples (3x3x3 mm<sup>3</sup>) are inserted into the oven via an object slide. Samples can either be illuminated from behind to obtain a silhouetted image, which is very effective for observing deformation processes and determining the contact angle, or they can be illuminated from above at the front - in this case, the surfaces are clearly visible and sintering processes can be followed.

The microscope has a magnification range between 3.5 and 11 and heating processes can be filmed with the attached CCD camera.

Filters can be added to enable very bright, glowing hot samples to be observed.

The oven and the camera are controlled via a computer program, which enables almost any heating ramp. Even periodical temperature profiles can be programmed.

**Fig. 2:** Programme for controlling the heating microscope.

### Specifications:

#### Measurable variables:

Different characteristic temperatures

#### Temperature:

Up to 1500 °C

#### Atmosphere:

Air

#### Magnification:

3.5x, 5x, 11x

#### Contact:

S. Vidi, Dipl. Phys.

Tel.: ++49-931/70564-50

Fax: ++49-931/70564-60

e-mail: vidi@zae.uni-wuerzburg.de

<http://www.zae-bayern.de>

#### Address:

ZAE Bayern

Am Hubland

97074 Würzburg

GERMANY