Emissionsgrad-Messanlage (EMMA)

Einleitung

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Emittance Measurement Apparatus (EMMA)

Introduction

Every body radiates due to its temperature, therefore emitting heat into the surrounding area.

Emittance is the measure of radiation emitted by a body due to its temperature compared to a black body. The value depends on the material, the state of the surface, the temperature T of the body, the wavelength λ and the direction θ in which radiation is emitted.

Due to the fact that emitted radiation increases significantly with increasing temperature, the emittance is important, for instance, for high temperature applications. Thermal insulation, solar collectors and thermal barrier coatings in turbines are examples where the emittance is a very important issue.

Measuring facilities

A black body and a sample oven are coupled to a FTIR spectrometer by way of an external additional optics tool. The **directional-spectral emittance** ϵ_λ is determined by measuring the spectral intensity emitted by the black body (reference radiator) and by the sample either vertical to the surface or at a particular angle according to:

$$\epsilon_\lambda(T) = \frac{i_{\lambda,\text{sample}}(T)}{i_{\lambda,\text{black body}}(T)}$$

Figure 1 shows the oven with a copper sample. The intensity emitted from the copper sample is detected by the spectrometer.


Fig. 1: The oven with a copper sample.

Heat-insulating layers in turbines

Partially yttria-stabilized zirconia (PYSZ), for instance, is used as a thermal barrier coating (TBC) in turbines.


Fig. 2: SEM image of a ZrO_2 layer deposited via electron-beam physical-vapour deposition (EP-PVD): On the left the surface, on the right a cross section.


Fig. 3: Apparent emittance of two layers of zirconia on steel compared to the emittance of the steel substrate alone.

Specifications

Measurable variable:	Emittance
Sample temperature:	200 °C ... 800 °C
Angle of radiation:	0° ... 70°
Wavelength:	2 µm ... 18 µm
Atmosphere:	vacuum, inert gas
Samples:	solid body, $\varnothing = 30 \dots 60$ mm

Contact:

Dr. J. Manara
Tel.: ++49-931/70564-46
Fax: ++49-931/70564-60
e-mail: manara@zae.uni-wuerzburg.de
<http://www.zae-bayern.de>

Address:

ZAE Bayern
Am Hubland
97074 Würzburg
GERMANY