# Hitzdraht-Methode

## Einleitung

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# Hot-Wire Method

## Introduction

Steady-state methods of measuring thermal conductivity often result in long measuring times. The advantages of the hot-wire method are the short measuring times and wide range of applications under the usual measuring conditions. Powders, solids, liquids and gases can be measured under the most diverse ambient conditions using the hot-wire method. Only relatively small samples are required.

## Measurement principle

A thin wire is embedded into the sample to be investigated, simultaneously serving as a heating element and a temperature sensor. During the experiment, the wire is heated with a constant electrical power source. The temporal development of the increase in temperature at the hot wire is calculated from the resistance of the wire. This temperature increase is primarily dependent on the thermal conductivity of the sample.

The thermal conductivity is determined by taking the thermal contact resistance between sample and wire as well as axial heat losses to the temporal temperature development into consideration.

## Application

ZAE Bayern has three hot-wire apparatuses designed for varying demands, which in essence differ in their range of operating temperatures. The total temperature range of all three apparatuses is -40 °C to 1600 °C, partially under a protective gas atmosphere. The gas pressure and external load can also be varied to a certain degree. Investigations to date have served, for instance, to optimize high-temperature insulation materials for use in space. The hot-wire procedure has also proved particularly suitable to thermally characterize powders used, for example, in vacuum insulation panels.

## Specifications

Measurable variable: thermal conductivity in the range between  $10^{-4} \text{ Wm}^{-1}\text{K}^{-1}$  and  $2 \text{ Wm}^{-1}\text{K}^{-1}$

Temperature range: -40 °C to 1600 °C

Atmosphere:  $10^{-6}$  hPa to 1000 hPa, different gases (nitrogen, argon, helium, air)

External load: 0 MPa to 0.1 MPa

Relative humidity: 10 % to 98 % within the temperature range between 10 °C and 95 °C

Samples: powders, liquids, gases, solids (approx. 100 mm long, diameter approx. 40 mm, or equivalent rectangular cross-sections)

### Contact:

F. Hemberger

Tel.: ++49-931/70564-26

Fax: ++49-931/70564-60

e-mail: [hemberger@zae.uni-wuerzburg.de](mailto:hemberger@zae.uni-wuerzburg.de)

<http://www.zae-bayern.de>

### Address:

ZAE Bayern

Am Hubland

97074 Würzburg

Germany