# Laserflash-Verfahren

## Einleitung

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# Laser-Flash Method

## Introduction


Transient heat transfer in solids is primarily determined by thermal diffusivity  $a$ , which can be directly derived using the laser-flash method. Thermal conductivity, the most important parameter in steady-state heat transfer, can also be determined using the laser-flash method. The thermal conductivity of a material can be calculated according to  $\lambda = a \rho c_p$  if the heat capacity  $c_p$  and density  $\rho$  are known.

Using the laser-flash method, thermal diffusivity can be precisely determined in a short time and even at high temperatures. A further advantage is that samples do not have to be specially prepared due to the non-contact connection of both the heating source and the temperature sensors.

A further significant application of the laser-flash method is determining thermal contact resistance on the boundary layers of composite materials, such as metal-to-metal bonds.

## Principle

With the laser-flash method, a laser beam serves as the heating source (cf. Fig. 1). The front of the disc-shaped sample under investigation is heated with a short laser pulse. The heat introduced in this way spreads over the whole sample also heating the back, whereby the time needed for this process is determined by the thermal diffusivity  $a$ . The heating process is recorded by an infrared detector. The time-temperature curve is mathematically analysed using numeric models to ascertain the thermal diffusivity  $a$ . It is possible to determine the temperature-dependent thermal conductivity by combining specific heat capacity measurements (using differential scanning calorimetry) and measuring the density of the sample.

**Figure 2:** Diagram of the laser-flash method

## Applications

The apparatus has two furnaces which enable measurements between room temperature and 1500 °C. Measurements can be carried out in air, inert gas or a vacuum. Examples of materials investigated to date include silicon carbide mono-crystals, aluminium nitride, ceramic coatings, glass, steel and plastic.

## Specifications

Measurable variables: Thermal diffusivity ( $0.5 \text{ mm}^2 \text{ s}^{-1}$  to  $200 \text{ mm}^2 \text{ s}^{-1}$ ), thermal contact resistance

Temperature: -10 °C to 1500 °C

Atmosphere: 0.01 to 1000 mbar, with different types of gas

Samples: Disc-shaped solids, thickness approx. 1 - 3 mm, diameter: 5 mm and 12.5 mm (for high-temperature measurements) or 5 - 24 mm in the range between room temperature and 250 °C.

### Contact:

F. Hemberger  
Tel.: ++49-931/70564-26  
Fax: ++49-931/70564-60  
e-mail: [hemberger@zae.uni-wuerzburg.de](mailto:hemberger@zae.uni-wuerzburg.de)  
<http://www.zae-bayern.de>

### Address:

ZAE Bayern  
Am Hubland  
97074 Würzburg  
Germany